

# SaliPen

THE HIGH-TECH MOUTH MOISTENER

## THE ONLY NATURAL THERAPY FOR DRY MOUTH!



**SaliPen** is a state of the art medical device that treats dry mouth naturally.



**SaliPen** is the only long-term dry mouth cure that makes your mouth wet **with your own saliva without any side effects or chemicals!**




**SaliPen** uses low-voltage electrostimulation to increase saliva production and is not felt by the patient. Compared with other dry mouth treatments - medications, gums, sprays, lozenges and others - SaliPen is also **very cost effective.**

**SaliPen** is worn typically for only a few minutes every day. With daily use, its effect lasts for hours. →

**SaliPen** consists of two flexible arms with electrodes that are placed in the mouth next to the inner gums of the lower wisdom teeth, underneath the tongue. The rest of the device stays outside the mouth.



A worldwide trial in 14 clinical centers showed significant improvement in oral wetness. SaliPen is approved by FDA (USA) and TGA (Australia) and the manufacturer, Saliwell Ltd., is ISO 13485:2016 certified.



### How to buy?

- 1 USA residents must obtain a prescription from a medical doctor or a dental practitioner and send it to [sales@saliwell.com](mailto:sales@saliwell.com)
- 2 Order at [www.saliwell.com/order-salipen-online](http://www.saliwell.com/order-salipen-online)